

HAUSSER SCIENTIFIC HEMACYTOMETERS

Hausser counting chambers are available in a variety of chamber types, for use in cell counting, examination of fluids and micro-organisms, and counting procedures on suspension emulsions. Hausser also offers accessories for use with counting chambers.

	CAT #:
SECTION I: DOUBLE NEUBAUER RULED METALLIZED HEMACYTOMETER	
Chamber only (bridges .1 mm deep)	3110
Chamber with two cover glasses (5000)	3100
Chamber with two cover glasses, pipettes and case (5000)	3120
Phase Counting Chamber with two cover glass (1461)	3200
Chamber with two cover glasses (5000)(no grid) (surface underneath is flat not concavity (hollow))	3102
SECTION II: ETCHED - TRANSPARENT FIELD WITH OPAQUE PATTERN LEVY DOUBLE NEUBAUER RULED HEMACYTOMETER	
Chamber only	3510
Chamber with two cover glasses (5000)	3500
Chamber with two cover glasses, pipettes and case (5000)	3520
OTHER HEMACYTOMETERS	
Palmer Counting Chamber (Nannoplankton) w/2 cg (square)	3850
McMaster Egg Slide (2 Cell chamber)	3875
Speirs-Levy Eosinophil Counting Chamber w/2 cg (5020)	3556
Fuchs-Rosenthal pattern - for cerebro-spinal fluid w/2 5000	3720
Sedgewick-Rafter Chamber w/1 cover glass (5040) no grid	3800
Howard Mold Counter with two cover glass (5080 & 5090)	3820
Petroff-Hausser Bacteria Counter w/2 cg (5051)(1 grid)	3900
Petroff-Hausser Bacteria Counter in case w/2cg (5051) (set)	3920
Petroff-Hausser Bacteria Counter (no grid) w/case & cg	3900NG
Petroff Hausser Bacteria Counter (10 microns) w/grid & 1 cg	3901
Petroff Hausser Bacteria Counter (20 microns)w/2grid & cg	3902
Petroff Hausser Bacteria Counter w/o grid (10 microns)w/cg	3951
Petroff Hausser Bacteria Counter w/o grid (20 microns)w/cg	3952
Nageotte, Bright Line, Double, Counting Chamber with cg (5400 & 5411)with cross flute for 2 different samples	4000
Nageotte, Bright Line, Double, Counting Chamber with cg (5400 & 5411) no cross flute	4001
SECTION III: ACCESSORIES FOR COUNTING CHAMBERS	
Cover glass 20 mm x 26 mm x .5mm	5000
Cover glass 25 mm x 28 mm x .5mm	5010
Cover glass 25 mm x 33 mm x .5mm (used on 3556)	5020
Cover glass 27 mm x 37 mm x .5mm	5030
Sedgewick Rafter cover glass 25 mm x 60mm x .5mm	5040
Petroff-Hausser reinforced cover glass 20mm x 26mm x .2mm	5050
Petroff-Hausser cover glass 20mm x 26 mm x .2mm	5051
Palmer thin cover glass (square)	5052
Whipple eyepiece micrometer w/grid used on 3800	5060
Howard eyepiece micrometer w/grid (21 mm diameter)	5070
Graduated Howard Mold Cover Glass, unlaminated	5075
Graduated Howard Mold Cover Glass, laminated	5075L
Howard thin cover glass 28 mm x 33 mm x .5mm	5080
Howard thick cover glass 28 mm x 33 mm x 1.0mm	5090
Nageotte cover glass 30 mm x 30 mm x .5mm (reusable)	5400
Nageotte cover glass 30 mm x 30 mm x 1.5mm (approx 35/box) (disposable)	5411
SECTION IV: LEICA (REICHERT) BRIGHTLINE HEMACYTOMETERS & ACCESSORIES	
Red Pipette	1456
White Pitette	1458
Pipette Filler	1450
Hemocytometer cover glass	1461
Phase counting chamber with two cover glass (1461)	1475
Bright-line counting chamber only	1490
Bright-line counting chamber with two cover glass (1461)	1492
Bright-line outfit complete	1483
Case only	1505